



In Technology,  
Simplicity is  
the Ultimate  
Sophistication

Introducing  
**EasyPro™ 3**  
Urodynamics



Urodynamic evaluations are fast becoming routine in the office environment, with only complicated cases referred to specialty centers. What accounts for this change in practice?

1. Third-party payers increasingly demand medical justification of surgical procedures - and for the treatment of patients with lower urinary tract symptoms (LUTS), that requires urodynamics.
2. Urodynamics is no longer complicated or cumbersome.

The EasyPro 3 is in the vanguard of a new generation in urodynamics technology, with a focus on ease-of-use and targeted clinical data that are easily interpreted. EasyPro is as easy to use as its name implies - yet offers features that no other system can match, like its unique electronic leak-point (ELP) detector.

The ELP detector is a unique, patented device that replaces hand-switches and event markers. It uses electronic sensors to accurately and automatically measure leak-point pressure. The ELP detector slips around any 6-8 Fr catheter and slides up to the meatus; at the first drop of urine, it reliably marks the event.

## Equipment Features

- ◆ 10-channel Urodynamics System: 3ch pressure, 2ch subtraction, 1ch electronic leak-point measurement, 2ch volume, 1ch sEMG and 1ch Uroflow
- ◆ Touch-screen console with integral PC and instrument
- ◆ Custom designed rolling stand for in-office portability
- ◆ Full-function remote controls pump, marks events, etc.
- ◆ Water-resistant, positive-connect cable system
- ◆ System complete with printer, pump, uroflow scale, commode and all necessary cabling

## Software Features

- ◆ Simple to use
- ◆ Automatic patient ID numbering, HIPAA compliant
- ◆ Infusion control from 1-100 mls
- ◆ Automatic catheter zero function
- ◆ Automatic and user defined parameters
- ◆ Summary study reporting (color) of relevant clinical data:
  - ◆ Area or line charting of all channels
  - ◆ All events marked with assoc. pressures/volumes/times
  - ◆ CMG/UPP/VLPP/Pressure-Flow/Uroflow reports including volumes, times, infusion rates, positions, etc.
  - ◆ Automatic calculation of residual bladder fluid
  - ◆ Nomogram, narrative, graphical and summary chart formats
- ◆ Store and retrieve functions for review/reprint of studies:
  - ◆ Add user notes and events
  - ◆ Save in standard electronic format
- ◆ Remote diagnostics capability
- ◆ On-line supplies ordering

## Technical Specifications

- ◆ Pressure ranges from 0-200 cm H<sub>2</sub>O
- ◆ EMG – Amplitude Range: 3-75  $\mu$ V, Frequency Range: 25-1,500 Hz, Input Impedance: 10 Meg $\Omega$
- ◆ Uroflow Range – 1.5 liter capacity

- ◆ 10-Channel Urodynamics System
- ◆ Patented ELP detector eliminates the guesswork, standardizes measurement of leak-point pressure
- ◆ Touch-screen simplifies operation
- ◆ Custom rolling stand for in-office portability
- ◆ Automatic reports for CMG, UPP, VLPP, Pressure-Flow, Uroflow, etc.



Uroflow



Remote-control

## Catheters and Supplies

SRS Medical offers a complete line of disposable catheter products and urodynamic supplies that features the latest designs and technology for adult, pediatric and specialty catheters, including water-charged, air-charged and solid-state (microtip) catheters.

## Related Products for Continence Care

SRS Medical is a worldwide leader in innovative medical devices for the diagnosis and conservative treatment of incontinence, a condition that affects the quality-of-life of over 30,000,000 adults. In addition to urodynamics, SRS Medical offers stand-alone uroflow, clinical biofeedback, neurostimulation and more. SRS Medical patient devices include StepFree™ vaginal weights, C3® penile clamps and the In-Flow™ female catheter (investigational device - not available in U.S.).

© 2003-2006. SRS Medical Corp. All rights reserved. P/N 21000345 Rev C

Headline paraphrased from Henry Wadsworth Longfellow (1807-1882). Original quotation is as follows: "In character, in manner, in style, in all things, the supreme excellence is simplicity."

Innovative Solutions for Incontinence and C3 are registered trademarks and EasyPro and StepFree are trademarks of SRS Medical Corp. In-Flow is a trademark of Personal Medical Corp.

**SRS**  
MEDICAL

Innovative Solutions for Incontinence®

800 345 5642 ◆ srsmedical.com