



## Eris™ UroScope™

# Bulking Agent Injector & Female Endoscope

- ◆ Patented rotatable sheath simplifies transurethral injections
- ◆ Can be used with all FDA-registered bulking agents
- ◆ Available in 0, 30 and 70 degree viewing angles
- ◆ Adaptable to ACMI, Storz and Wolff fittings

# Eris UroScope

## Bulking Agent Injector & Female Endoscope

The Eris UroScope can be used with all FDA-registered bulking agents. Its new rotatable Needle-Guide Sheath (patented), for exclusive use on the 0 degree Female Endoscope, greatly simplifies transurethral injection procedures.

Recent data concerning the use of bulking agents for the treatment of female SUI indicates that the best long-term results are obtained via a circumferential injection pattern. The Eris UroScope is the best tool for accomplishing this. Unlike other injectors that require manipulation and support of the cable each time you change sites, the Eris UroScope has a rotatable sheath, allowing you to change position easily and reliably.

The Eris UroScope is also an excellent female endoscope, designed by Gynecologists for ease of use and patient comfort. It is available in 0, 30 and 70 degree viewing angles (Urethroscopy, Hysteroscopy and Cystoscopy) and is adaptable to ACMI, Storz and Wolff fittings.

### *Related Products for Continence Care*

SRS Medical is a worldwide leader in innovative medical devices for the diagnosis and conservative treatment of incontinence, a condition that affects the quality-of-life of over 30,000,000 adults. In addition to the UroScope, SRS Medical offers urodynamics, clinical biofeedback, neurostimulation and the In-Flow™ device, an advanced technology for women with atonic bladder. (Not available in U.S.)

© 2003 SRS Medical Systems, Inc. All rights reserved.

U.S. Patent No. 6,053,860. Innovative Solutions for Incontinence is a registered trademark and Eris, UroScope and In-Flow are trademarks of SRS Medical Systems, Inc.



*Innovative Solutions for Incontinence®*

**800 345 5642 ♦ [srsmedical.com](http://srsmedical.com)**