

In-Office  
Urodynamics  
Made Easy

EasyPro™ 4



# EasyPro™ 4 Urodynamics



With the Only Electronic Leak-Point Detector

- 12-Channel Urodynamics System
- Fast, accurate measurement of leak-point using electronic detector (an SRS Medical exclusive)
- Easy-to-use touch-screen interface
- Custom rolling stand for in-office portability
- Informative, easy-to-understand reports

Urodynamic evaluations are fast becoming routine in the office environment, with only complicated cases referred to specialty centers. What accounts for this change in practice?

1. Third-party payers increasingly demand medical justification of surgical procedures – and for patients with lower urinary tract symptoms (LUTS), that requires urodynamics.
2. Urodynamics is no longer complicated or cumbersome.

The EasyPro4 is in the vanguard of a new generation in urodynamics technology, with a focus on ease-of-use and targeted clinical data that are easily interpreted. EasyPro4 is as easy to use as its name implies – yet offers features that no other system can match, like its unique leak-point detector.

## Disposable Electronic Leak-Point Detector

The point at which a patient first leaks urine is one of the most clinically significant measurements in urodynamics. Historically, these measurements have depended on the clinician's visual acuity and reflexes. Now you can take the uncertainty out of leak-point. The patented LPD™3d leak-point detector (top right) allows you to make reliable measurements every time.

The LPD-3d is a unique, patented device that replaces hand-switches and event markers. It uses an electronic sensor to accurately and automatically measure leaks. The LPD-3d slips around any 5-10 Fr catheter and slides up to the meatus; at the first drop of urine, it reliably marks the event. In addition, only the LPD-3d can reliably detect repeated leaks.

## Advanced Urodynamics Made Simple

The EasyPro4 features a simple, consistent touch-screen interface that makes it easy to operate. At the end of a urodynamics study, the system prints out a comprehensive, but easy-to-understand set of reports that includes data for CMG, VPS, UPP, uroflow, etc.

## Everything You Need Now

The EasyPro4 includes everything you need for in-office urodynamics:

- ◆ EasyPro4 instrument console with 4ch pressure, 2ch EMG, direct-connect for leak-point detector;
- ◆ Rolling 2" pole cart, 17" touch-screen PC, printer, and hospital-grade power conditioner;
- ◆ EasyFlo™4 wireless uroflow scale and commode;
- ◆ Peristaltic pump with primer control;
- ◆ Full-function remote control with all event markers, pump control, etc.; and
- ◆ All necessary cabling, including for T-Doc brand air-charged catheters



Wireless Uroflow



Air-Charged Catheters

## And Everything You'll Want Later

The EasyPro4 incorporates two channels of EMG, allowing use for pelvic floor rehabilitation with optional iMP™ software. In addition, as part of the Clarity™ uDx urodiagnostic workstation, the EasyPro4 has a built-in upgrade path. Options include:

- ◆ UroScan® Ultrasound for abdominal imaging and bladder volume measurements,
- ◆ The CT3000 UroCuff™ Test for non-invasive male voiding pressure, and
- ◆ Additional EasyFlo4 wireless uroflow scales.

## Catheters and Supplies

SRS Medical offers the most competitive pricing on urodynamic supplies, including pump tubing, EMG electrodes and all types of catheters. Call us.

© 2003-2009. SRS Medical Corp. All rights reserved. P/N 21000345 Rev E

UroScan is a registered trademark and EasyPro, LPD, EasyFlo, iMP, Clarity and UroCuff are trademarks of SRS Medical Corp.



Pelvic Health Technologies

800 345 5642 • srsmedical.com