



Your Complete Source for Electrodes and Sensors*



SenseRx™ Vaginal and Anal Sensors

Product Nos. 4509C (Vaginal) and 4510C (Anal)

SenseRx sensors are the only internal sensors with extensive documentation of clinical efficacy and have been cited in numerous medical texts and peer-reviewed journal articles. With the best “fit & finish” of any internal sensor, these remain the Quality Standard. Internal sensors make the most faithful measurements of EMG, since they directly contact the pelvic floor muscles. Their longitudinal (lengthwise, not circular) detection electrodes maintain the proper orientation to muscle fibers for highest accuracy. Because they are self-seating, SenseRx sensors can be patient-inserted and patients can remain fully clothed during treatments. Can be used for both EMG and pelvic floor stimulation.



Eisman-Tries MEP™ Dual Anal Sensor

Product No. 4520

The MEP monitors the EMG activity of the levator ani muscle, including simultaneous measurements of the proximal and distal bundles of the external anal sphincter (EAS). The deep measurement of the EAS makes the MEP uniquely suited for post-prostatectomy pelvic muscle rehab, since it allows training to be specifically directed to improving tonic support of the bladder neck.



IVS-2™ Intravaginal Sensor

Product No. 4525

The IVS-2 is the smallest, most comfortable female EMG sensor available and features gold detection electrodes. Gold is a non-polarizing material that provides consistently superior measurements and is also hypoallergenic. The IVS-2 is inserted similarly to a tampon and positions above the urogenital diaphragm. Its small size and unique placement allows it to measure muscle activity even when the patient is moving. This is an important benefit, since certain movements are closely associated with urine loss.



Surface EMG Electrodes

Product No. 4630 (Sets of 3)

SRS Medical surface EMG electrodes provide excellent signal stability and are optimized for patient comfort and safety. They are soft, flexible, latex-free and their hydrogel adhesive is hypo-allergenic.



EMG Electrodes with Wires

Product No. 4639

Same as the SRS Medical surface EMG electrodes, but with the addition of integral wires, which makes it easier to meet infection control standards for perianal placement. The entire electrode assembly is disposable - you never have to clean cables again.

*All electrodes and sensors can be used with most commercially available instruments. Adaptor may be required for non-SRS Medical instruments.

© 2000-2010. SRS Medical Corp. All rights reserved. SenseRx, IVS-2 and MEP are trademarks of SRS Medical Corp. ESDS-01 Rev. B



Advancing Pelvic Health

800 345 5642 • srsmedical.com