

The Clarity™ uDx

Urodiagnostic Workstation

Now...
One System
Does It All



— Touch-Screen Operation
A simple, consistent interface
for all diagnostic tests

— Clarity Uro-Software
One database provides single
EMR node for all test data

UroScan® Ultrasound
Bladder Volume and Imaging

CT3000 UroCuff™ Test
Non-Invasive Voiding Pressure

— EasyPro™ 4 Urodynamics
with optional iMP™ EMG Rehab

Turnkey System
Includes stand, printer,
power isolator, etc.

EasyFlo™ 4 Uroflow
Wireless Scale/s




6 Urodiagnostics

1. Urodynamics
2. Non-invasive male voiding pressure
3. Ultrasound for abdominal imaging
4. Uroflow
5. Ultrasound for bladder volume
6. EMG for Dx and pelvic floor rehabilitation

1 Workstation

- One simple, consistent interface
- One database
- One EMR node
- One compact instrument cart

Choose Clarity™ uDx (on front) or
Start with a Dedicated System and Build Up:




CT3000 UroCuff™ Test

Includes:

- Clarity base station with 2" pole cart, 17" touch-screen PC, printer, hospital-grade power conditioner
- CT3000 instrument console
- EasyFlo™4 uroflow, commode

Options:

- UroScan® ultrasound
- CT3000 2ch EMG
- Additional EasyFlo4 uroflows




EasyPro™ 4 Urodynamics

Includes:

- Clarity base station with 2" pole cart, 17" touch-screen PC, printer, hospital-grade power conditioner
- EasyPro4 instrument console
- EasyFlo4 uroflow, commode, pump

Options:

- UroScan ultrasound
- iMP™ software for pelvic floor rehab
- Additional EasyFlo4 uroflows



VETS™ Non-In Dx/Tx

Includes:

- Clarity base station with 2" pole cart, 17" touch-screen PC, printer, hospital-grade power conditioner
- Aware™ EMG dual muscle monitor
- EasyFlo4 uroflow, commode
- UroScan ultrasound

Options:

- Additional EasyFlo4 uroflows